


**IPL HAIR REMOVAL HANDSET
FOR HOME USE**

IPL-1



User Guide

Thank you for choosing our IPL handset. To achieve the ideal results, please read through this user guide carefully before using your IPL handset.

1. Safety Tips

To protect your skin, this handset is not suitable for the following conditions:



- You have sensory disturbance of skin.
- You have skin diseases, such as allergic dermatitis, eczema etc
- You are affected by photo-sensitization.
- You are menstruating, pregnant or lactating.
- You have sunburn.
- Broken skin; a scab or healed wound; skin cancer or hemangioma.
- Do not use on or near your eyes or any artificial aesthetic, prosthetic or plastic surgery area.
- Do not use the handset if you have dark spots such as large freckles, moles or blisters on the treatment area.

2. Product Information

1) What's in the box:

1. IPL Laser Hair Removal Handset
2. Power Adaptor
3. User Guide
4. Warranty Card

2) Handset Description



3. Technical Parameters

Lamp Tube	Xenon Quartz Tube
Technology	Intense Pulsed Light
Energy Level	1.5- 3.9 J/ cm ²
Focus Size	33*10mm ²
Wavelength	470nm ~1100nm
Adaptor	100~240V AC Input 12V3A/2.5A DC Output
Operating Temperature	5~30 °C
Weight	186g
Size	124*83*48.5mm ³

4. Functions

1) Hair Removal

IPL Hair Removal Handset can be used on any body part : Face; Arms; Legs; Underarms; Brazilian Area. Please check the Skin Tone and Hair Tone Tables below to determine whether IPL is suitable for you.

Skin Tone

SKin Type	Description	Recommended Intensity Level
	Ivory white skin that burns easily and never tans. Usually very light skin tone and light eyes	5
	White skin that burns easily and tans minimally with difficulty.	5
	White skin that burns moderately and tans moderately and uniformly.	4
	Beige-olive,lightly tanned skin that burns minimally and tans moderately and easily.	3~4
	Moderate brown or tanned skin that rarely burns and tans profusely	1~2
	Dark brown or black skin that never burns or tans profusely. Usually very dark skin tone, and dark eyes.	Not suitable

Hair Tone

	White/Grey	no
	Light Gold	no
	Red	yes
	Light Brown	yes
	Brown	yes
	Dark Brown	yes
	Black	yes

5. Instructions

- 1) Clean the flash window with the alcohol wipe.
- 2) Plug the handset in and press the green button on the back to turn it on.
- 3) The green button is also used to switch through the intensity levels (1 for low - 5 for high).
- 4) The blue lights on the front will indicate the intensity level.

- 5) Complete your spot testing before using the handset.
- 6) For the best results, you should use your handset on the highest possible intensity level that is not painful. Gradually increasing the intensity level after each session if the experience was not painful.
- 7) The handset has a sensor, and will not flash unless the flash window is in direct contact with your skin. If is not in direct contact with your skin, the red light will stay illuminated and the handset will not flash.
- 8) To use, ensure that the light window is completely flat against the skin.
- 9) Indicator light screen will begin to blink once you' ve placed the light window onto your skin. This indicates that the handset is ready to be used.
- 10) Press the flash button.
- 11) For stubborn areas, you can treat the area with multiple passes (running your handset over the same area up to 2 times).
- 12) Press and hold the green button on the back to turn it off.
- 13) To clean the handset, dampen a cloth with a few drops of water and wipe over the light window. Make sure that the IPL is unplugged from the wall.
- 14) If the handset inactive for 5 minutes, the handset turns into sleep mode, gently press the power button to wake it up. When restart, the intensity level back to level 1 automatically.

Note:

After frequent use of the handset, there will be few small white spots inside flash window. This is normal and will not affect the use and lifespan of handset.

6. Auto mode/Manual mode switching

The default mode is manual, press the flash button for 2 seconds to switch handset to auto mode. Under Auto mode, the handset automatically flashes when detected the skin. The flash stops when handset away from skin. To quit auto mode, please long press the flash button for 2 seconds, then the handset back to manual mode.

7. Overheat protection

The overheat may occurred if the handset continuously use for a long time. In this situation, and the handset cannot flash. When the overheat protection starts, the radiator fan is still working inside to help cool down the temperature. Please allow your handset take 2-5 minutes to cool down before continuing use.

After the internal temperature returns to normal, the handset starts over, but at this time the intensity level switches to level 1, if you need the higher intensity level, you can reset it by yourself.

8. Spot Test Before Use

We recommend you complete a spot test before starting your full body treatment.

The spot test is intended to determine your skin's suitability to IPL, and to help you select the optimal energy level for each treat area.

- 1) Do a spot test on each area you intend to treat 48 hours prior to the full treatment.
- 2) Prior to use, make sure your skin is shaved, clean and dry (free of any residues such as powders, creams, make-up, antiperspirants or deodorants).
- 3) Set up the device following the instructions above.
- 4) Ensure flash window is completely flat against your skin.
- 5) Deliver the first flash on in intensity level 1 by pressing the flash button.
- 6) Remove the device from skin and place it on the next spot of the test area.
- 7) Use power button to increase the intensity level by 1. The power lights turn up according to the selected energy level. Deliver the second flash on intensity level 2 by pressing the flash button.
- 8) When you feel comfortable, please repeat steps 5-7 until you have reached the maximum intensity level recommended for your skin tone. Stop the test if you notice any adverse reactions.

- 9) Wait 48 hours and examine the tested areas. If the skin in the test area seems normal, you can proceed with a full treatment for those areas.

9. Cautions

- 1) Do not aim the handset or flash window at your eyes.
- 2) Do not use the device around eyes, eyebrow or eyelashes as this could cause serious and permanent eye damage
- 3) It is normal if your skin gets redness after IPL treatment. You can cool the skin with ice water or cold wet towel.
- 4) Do not clean the flash window directly with alcohol spray.
- 5) Do not use with any Gel cosmetic products.
- 6) Try not to eat the photographic food, such as spinach, bok-choy, Chinese broccoli, etc.
- 7) For the best results, you can use the acceptable higher intensity level which not hurting yourself.
- 8) Avoid long exposure in direct sunlight (including tanning beds) for at least 2 weeks prior and 48 hours post treatment. If you are exposed to the sun after treatment, we advise using an SPF sunscreen to protect your skin.
- 9) After the IPL treatment, please moisturizing your skin with non-alcoholic cream or lotion.

- 10) Do not use in or around water.
- 11) For the darker skin tone users, the skin reactions may appear slowly.
Please watch the skin condition before increasing the intensity level.
- 12) Do not cover the air inlet and the air vents in IPL treatment.
- 13) The remained hair on the handset may causes cosmetically darken the handset, but it will not affect the normal use of the handset.
- 14) It is normal the lamp tube be darken in its late life that will not affect the using of the handset.
- 15) Avoid storing in places with strong magnetic fields and severe humidity and dust or high temperature places.
- 16) Keep out of reach of children. We recommend that the handset should not be used by anyone under the age of 16.
- 17) Do not open or attempt to repair the device at any time.
- 18) Do not strongly strike or knock the handset.
- 19) Do not use the handset for any purpose other than the designated function.